
# KRUUSE PVC Endotracheal Tubes

With cuff

Without cuff

## KRUUSE PVC Endotracheal Tubes

### Instructions for use

Read the entire instructions before using the product.

KRUUSE PVC Endotracheal Tubes are made of polyvinyl chloride.

Regarding the tubes with cuffs, the filling system consists of a filling tube, a pilot balloon, a self-closing valve with a Luer connection. All KRUUSE endotracheal tubes are with Murphy's eye and in a curved, anatomically shaped form. The cuffed endotracheal tubes are HVLP-type (High Volume, Low Pressure).

### Instructions for intubation

1. Carefully remove the sterile endotracheal tube from the pouch
2. The passage of the main tube lumen should be checked prior to intubation
3. If the tube has to be cut (shortened) prior to intubation, remove the fitted connector and cut the tube at the appropriate point. The connector should then firmly be reinserted into the endotracheal tube

### The following only regards endotracheal tubes with cuffs

4. Check the cuff, pilot balloon and valve on its leak-proofness and integrity prior to use. Insert a Luer syringe into the Luer adapter for the filling system and inject enough air to fully inflate the cuff
5. After the test-inflation of the cuff, completely evacuate the air prior to intubation
6. Intubate the patient in accordance with clinical practice noting the specific Precautions (cuff-related) below
7. Once the patient is intubated, use a Luer syringe to inflate the cuff with just enough air to provide an effective seal between the endotracheal tube and trachea. The use of minimal occluding volume, minimum leak techniques and continuous monitoring of the cuff pressure can help reduce the incidence of many of the adverse reactions associated with the use of cuffed tracheal tubes
8. Remove the syringe from the valve after cuff inflation. Leaving the syringe attached will keep the valve open, permitting the cuff to deflate
9. Fix the tube to the patient in accordance with clinical practice
10. Check that the inflation system is not leaking. Cuff pressure should be measured periodically or permanently
11. Prior to extubation, deflate the cuff by inserting a syringe into the cuff filling system and removing all of the air contained within the cuff until the pilot balloon has collapsed
12. Extubate the patient in accordance with clinical practice. Dispose the endotracheal tube

### Precautions (cuff-related)

- Inflation of the cuff by "feel" alone or by using a measured amount of air is not recommended, since resistance is an unreliable guide during inflation.
- The use of a cuff pressure gauge (AG Cuffil) will help to monitor and adjust the cuff pressure
- Do not overinflate the cuff. The cuff pressure should normally not exceed 25 cm H<sub>2</sub>O. Over-inflation can result in tracheal damage, rupture of the cuff with subsequent deflation, or cuff distortion, which may lead to airway blockage
- The use of lidocaine topical aerosols might result in formation of pinholes in the cuff
- Use of a water soluble lubricant gel is recommended

- Various bony anatomical structures (e.g. teeth) within the intubation routes or any intubation aid with sharp surfaces can damage the integrity of the cuff. Care must be taken to avoid damaging the thin-walled cuff during intubation. If the cuff is damaged, the tube should not be used.
- Deflate the cuff prior to repositioning the tube. Movement of the tube with the cuff inflated could result in damage to the cuff or in patient injury

### General Precautions

- In case the tube is being cut to length, when reinserting the connector, ensure that it is not reinserted in the proximity of the joint between the inflation system and the tube
- The user should be aware of anatomical variations including the length of the airway
- When the position of the tube or of the patient is altered after intubation, it is essential to verify that the tube remains correctly positioned
- If the tracheal tube is lubricated prior to intubation, it is essential to verify that the lubricant does not enter and occlude the tube lumen or the cuff filling system, thereby preventing ventilation, or causing damage to the cuff
- Use of lubricating solutions to ease reinsertion of the connector is not recommended as it may contribute to accidental disconnection
- The device should be disposed of in accordance with applicable national regulations

### Single use only

Re-use or re-sterilization of devices intended for single use only may result in degraded performance or loss of functionality and/or in exposure to viral, bacterial, fungal, or prionic pathogens.

### Sterility

Sterile unless packaging has been opened or damaged. Discard after single use, do not re-use nor re-sterilize. Do not use if the product packaging is compromised.

### Storage

Exposure to elevated temperature and UV light and excessive moisture should be avoided during storage. Keep away from sunlight and keep dry.

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